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CLSYRB

Seayu is a registered trademark of Swiss Eyewear Group (International) AG, Freilagerstrasse 39, 8047 Zürich, Switzerland

EU: Swiss Eyewear Group (Italia) S.r.l., Via Piave 11, 21051 Arcisate (Varese), Italy

Product complying with EN 14139:2010 and Reg. EU 2017/745

More information and the Declaration of Conformity are available on iseayu.com

seayu.

READING GLASSES
PREMIUM BLUE LIGHT FILTER

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glasses is 61 mm. These reading glasses are suitable for reading and near vision only and are not recommended for distance vision and are not to be used for eye protection. Solely an examination by a qualified optical professional can fully determine the precise prescription for optimal vision. An optical professional can further check on the general health of your eyes. These reading glasses are not to be used when driving a motor vehicle. Clean glasses with water and mild soap. Dry with a soft cloth. Avoid scratching the glasses. Keep in case. Replace the glasses if there is noticeable damage. **INFORMATION REGARDING BLUE LIGHT** Electronic devices such as smart phones and computer screens emit blue light. Exposure to blue light can cause eye strain or sleeplessness. Recent studies have also linked exposure to high energy visible light (HEV) to more

serious eye diseases. This bluelight filtering technology significantly reduces bluelight in the HEV spectrum of 380 to 430 Nanometers. **WARNINGS, PRECAUTIONS AND SIDE EFFECTS** These bluelight filtering glasses do not protect eyes from sunlight and should also not be used by children and infants. The use of these bluelight filtering glasses may lead to changes in the sleeping patterns, such as falling asleep easier. If the glasses have a negative effect on sleep patterns or any other adverse effects, immediately stop using the glasses. In low lighting conditions the use of these bluelight filtering glasses may reduce visual performance. Eye strain may be caused by a range of factors beyond exposure to the bluelight emitted from electronic devices. Causes may include refractive errors and the distance of the screen from the eyes, as well as conditions in the

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