


ORDO



3 Floss Tips
Orthodontic, Periodontal
& Premium Clean
Compatible with:
Ordo Hydro Sonic
Water Flosser



Orthodontic Tip
Tapered bristles help clean around braces,
bridges, implants and other dental work.



Periodontal Tip
Designed to rinse periodontal pockets,
helping to protect against gum disease.



Premium Clean Tip
Ideal for excellent everyday cleaning.
Suitable for most people.



We recommend replacing
your floss tips every 6
months.

Replace sooner if tips
become worn or dirty.

Directions

- Remove the old floss tip by pushing the release button on the back of your water flosser.
- Slide the new floss tip into place until it clicks.
- Rinse your floss tip before first use and rinse thoroughly after each use.

For more information on Ordo's products
and services please visit ordolife.com

Authorised Rep Compliance Ltd., Ground Floor,
71 Lower Baggot Street, Dublin, D02 P593, IE.
Ordo, Manchester, WA14 4DR, UK.
Made in China



Batch: 2310-SG3



ORDO



Tough on plaque
gentle on gums



Orthodontic Tip
for dental work



Periodontal Tip
for gum health



Premium Clean Tip
for a deep clean



3 Pack

HYDRO SONIC
FLOSS TIPS

Dental Selection

ORDO