hairmax®
TRANSFORMATIVE HAIR REGROWTH

USER MANUAL

ULTIMA **12**
LASERCOMB

Hair Growth **Laser Light Device**

INTRODUCTION

The Hairmax® laser device is a home-use laser device that utilizes clinically proven technology to treat male and female pattern hair loss, also known as hereditary/hormonal hair loss in men and women. The Hairmax® device is an over-the-counter (OTC) device that can be sold directly to the public.

Indication

Hairmax® laser devices are indicated to promote hair growth in males with Androgenetic alopecia who have Norwood-Hamilton Classifications of IIa - V, and in females with Androgenetic alopecia who have Ludwig (Savin) Classifications I - II or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I - IV.

How Hairmax® laser devices work

Hairmax® laser devices work by providing therapeutic laser energy to stimulate hair follicles. Hairmax® laser devices should be used 3 times per week, spaced out every other day. For example, Monday, Wednesday, Friday. The Hairmax® laser device should not be used more than the recommended treatment time per day. Clinical studies demonstrate that most users see benefits in as little as 16 weeks.

Hairmax® device results

Hairmax® devices provide laser therapy that has been proven safe and effective in 7 clinical studies with over 460 men and women. In these clinical studies, over 93% of the participants experienced significant hair growth, on average, over 20 new hairs per square centimeter/129 per square inch.

It is important that you read and understand the content of this user manual before using the Hairmax® laser device.


WHO IS THE BEST CANDIDATE

How do you know if the Hairmax® laser device is right for you?

Hairmax® laser devices can be used by men and women with thinning hair (including color treated hair) or pattern baldness caused by a hereditary/hormonal condition. Doctors use a system known as the Norwood Hamilton Classification (men) and the Ludwig-Savin Classification (women) to describe the degree of hair loss. Use the scales below to determine if Hairmax® laser devices are right for your hair loss. Hairmax® laser devices utilize clinically proven technology to treat hair loss as shown in the dark color areas below.

Classification for men and women

Norwood Hamilton Classification - Males

Ludwig-Savin Classification - Females

PARTS**WHAT IS INCLUDED**

1. Ultima 12 LaserComb
2. Charging cradle/stand
3. Power supply with cord
4. Quick Start Guide

Notice: Entire Hairmax® laser device assembly is considered as a Type BF applied part. The Hairmax® laser device is manufactured in an ISO 13485 certified facility using the highest quality components. Based on the durability of the lasers and electronic components, the expected life of the laser device is designed to last up to 10 years of regular use. The battery is rated for at least 300 cycles and may need to be replaced by manufacturer during the product expected life cycle. Please refer to the warranty section for additional details on warranty.

GETTING STARTED

Read this manual completely before you start to use your laser device. The instructions in this user manual are intended for individuals with a minimum of an eighth grade education. The device is safe for home use. The consumer is the intended operator of this equipment.

Use the Hairmax Ultima 12 LaserComb for 8 minutes per treatment, three treatments per week. The treatments should be spaced every other day, for example, Monday, Wednesday, Friday. Do not use the laser device more than the recommended treatment time. Most users see benefits in as little as 16 weeks (based on clinical studies).

When you receive the laser device, take all the parts out of the case. You must charge the device for 2-3 hours before initial use.

3 times
per week

TESTING FOR SENSITIVITY

Prior to treatment, test your skin for sensitivity by shining the laser light on the inside of your forearm for 2 minutes. Make sure your forearm is clean before testing. If you see any reaction, do not use the laser device and contact customer service.

CHARGING OPTIONS

OPTION 1: Charge the laser device using the charging stand and power cord.

Prior to charging the laser device, make sure the device is turned off and the teeth are facing away from you.

Plug the power cord into the charging cradle/stand and the other end into an electrical wall outlet.

The indicator light on the power button will flash during the charging process. Once charging is complete, the light will stop flashing.

OPTION 2: Charge the laser device by placing the power cord directly into the laser device and the other end into an electrical wall outlet.

The device indicator white light will flash during the charging process. Once charging is complete, the light will stop flashing.

DURING USE The indicator light will indicate the battery level. A steady, non-flashing light indicates the battery is full. As the battery power diminishes, the light will blink slowly. When the battery is getting low and needs to be recharged, the indicator light will start to flash rapidly. When this happens, charge the device as described above.

If the laser device needs to be recharged during your treatment, allow it to fully charge before continuing.



NOTE – Your Hairmax® laser device cannot be used while charging.

Before Starting Your Treatment:

1. Make sure the battery is **FULLY CHARGED**.
2. Test your Hairmax® laser device for sensitivity.


TREATMENT METHOD

Turn on your device by pressing the power button (lower button) and hold until the power button light turns on. Next, press the laser button (upper button) until it lights up. That means your device is ready for use.

After you have turned on the laser device, place it flat on your scalp so that both rows of the teeth are touching your scalp.

Start at your hairline. Move the laser device from front to back of scalp in a slow combing motion, pausing every ½ inch (1 cm) for 4 seconds. Do not lift the device off your scalp while combing from front to back. Next, lift the device off your scalp and start at a new position at the hairline.


Continue in this manner from front to back and then side to side until you cover your entire scalp with multiple passes. Repeat these passes for the total treatment of 8 minutes.

The information display on the top of the laser device will show the elapsed treatment time. The outer circle will represent the number of seconds of usage. The number on the middle represents the number of minutes that you have been using your device.

This example shows that you have been using the laser device for 4 minutes and 45 seconds.

If you have long hair, use your fingers to lift your hair out of the way while using the laser device to be sure that the laser energy reaches your scalp.

If you wish to turn off the beep/vibration during the treatment, press the top laser button once. The device and laser will remain on. To re-activate beep and vibration, simply press the laser button twice.

After treatment, turn off the laser device by pressing and holding the power button (lower button) until the indicator light turns off. Store the device in the charging cradle for future use. You may leave the power cord plugged into the charging cradle at all times.

The Ultima 12 LaserComb will turn off automatically after 15 minutes.

Please Note:

It is important to move the laser device very slowly. Remember, you are bathing the scalp with laser light. If you move it too fast, the treatment will not achieve the optimal results to promote hair growth.

TROUBLESHOOTING

If you are experiencing difficulties with your Hairmax® laser device, please review the tips below to resolve common issues. If further assistance is needed, please visit our Hairmax.com FAQ or contact Customer Support info@hairmax.com.

If the laser device will not turn on

- Ensure the device is fully charged.
- The device cannot be powered on while plugged in.

If the laser device turns off during use

- Confirm the device is fully charged.
- Press and hold the power button for 4 seconds, then press the laser button.

If the laser device will not charge

- Verify you are using a functioning electrical outlet.
- Ensure the power cord is inserted correctly into the device or charging stand.
- Always use the original charger that came with your Hairmax® laser device.

If the device vibrates irregularly or vibrates without laser light

- Confirm the device is fully charged.
- Ensure you are using the original charger provided with the unit.

If the laser device will not turn off

- Press and hold the power button for 4 seconds.

▲ **IMPORTANT:** Any attempt to self-serve the laser device, unless directed by Hairmax Customer Service, will void your warranty. For more information, troubleshooting support, or product guidance, please visit our Hairmax FAQ page or contact our customer care team directly.

LASER PARAMETERS

Wavelength: 655nm +/- 10nm visible red light

Laser power for classification: <5mW CW

Divergence: 57 mrad

Number of lasers: 12

Laser MTBF (mean time between failures):

10,000 hours at 25°C ambient temperature

This laser product is designated as Class 3R during all procedures of operation.

Note: Hairmax® laser devices have no user-serviceable parts.

USE CONDITIONS

Operating temperature/altitude/atmospheric pressure/humidity:

5°C to 35°C/0 to 2,000 m/80kPa to 106kPa/ 0% to 90%, relative humidity

Storage & transportation temperature/humidity/atmospheric pressure:

-20°C to 60°C / 0% to 90% / 80kPa to 106kPa

Ingress of water or particulate matter: IP22

Battery charging cycles: 300 cycles

ROHS (restriction of hazardous substances): 2011/65/EU

Product warranty: 1 year

Age of use: Adult

Device size and weight: 25.1 x 6.1 x 5.1 cm / 9.8 x 2.4 x 2.0 in / 184 grams / 0.4 lbs

Cradle size and weight: 11 x 6.6 x 14 cm / 4.3 x 2.6 x 5.5 in / 181 grams / 0.4 lbs

POWER SUPPLY SPECIFICATIONS

Output = 5.0 V DC \pm .25 V DC / 1.0A

The power supply for your Hairmax® laser device can be used in the United States and Internationally. This power supply will adapt to all power supplies from 100 volt/60 Hz (US) to 240 volt/50 Hz. If outlets in your country have different plugs than the power supply that comes with this laser device, use an international plug adapter (not supplied). Only use the power supply and power cord provided by the manufacturer, any other power supply or power cord may damage your device.

WARNINGS

- ▲ Do not stare directly at the red laser light or its reflection in a mirror because it can temporarily irritate your eyes. If you stare at the light for too long, it could harm your eyes. Also, never use a magnifying glass when using the laser device because it can cause temporary irritation to your eyes.
- ▲ Keep the laser device out of the reach of children, pets and pests.
- ▲ Do not drop the laser device in water because you may get an electrical shock. If the device gets damaged by water, contact customer service.
- ▲ If the cord is damaged on the laser device, do not use it because you could get an electrical shock. If the cord gets damaged, contact customer service.
- ▲ Do not use the laser device near any heated surfaces because this could cause the device not to work properly and may cause an electrical shock, contact customer service.
- ▲ Laser light - avoid direct eye exposure.
- ▲ This is a harmless low-level ESD (Electrostatic Discharge) sensitive device. Please do not touch the charging contacts or USB-port. Touching the contacts or USB-port is not of the essential performance and will not cause a basic safety issue. If touching the contacts or USB-port causes the device to malfunction, reset the device by powering off and back on.
- ▲ No modification of this equipment is allowed.
- ▲ If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use. Contact customer service.
- ▲ In case of a serious incident, seek medical attention immediately. Contact manufacturer and local authorities to report any serious incidents.
- ▲ If you are experiencing any allergic reaction from handling or using the laser device, stop using it immediately.
- ▲ Only use the power supply and power cord provided by the manufacturer, any other power supply or power cord may damage the laser device.

RISKS

The Hairmax® laser device has been clinically tested on both men and women. In the clinical study for men, 5% of the men said they had itching or tingling from using the Hairmax® laser device. None of the women in the clinical study reported any of these side effects. If you have itching, tingling or any other side effect that bothers you, stop using the device and contact customer service.

CAUTIONS

- ⚠ If you are taking any medications or products that make you sensitive to light, test the laser device to make sure you do not have a reaction. Some of the things that can cause increased sensitivity to the device could be medicines for colds, allergies, pain medicine and drugs to treat infections. To test your skin for sensitivity, turn the device on, hold the device on one spot, shining the light on your forearm. After two minutes, turn the device off and set aside. Look at the spot where you held the laser light over the next five minutes. If you see any reaction such as redness or rash, do not use the device and contact customer service for information on how to return the product and get a refund.
- ⚠ Only use the laser device as directed. Do not try to change the way the device operates. Altering controls, adjustments, performance or procedures other than those specified herein may result in hazardous radiation (laser) exposure. Contact customer service.
- ⚠ Never use acetone or other solvents on any part of the laser device. Use of solvents on the window can block the laser light from coming through.
- ⚠ It is suggested to charge the laser device every 3 months even if it is not used.
- ⚠ For assistance setting up, using or maintaining the laser device, or to report unexpected operations or events, contact customer service.
- ⚠ Do not place the power adapter where it is difficult to disconnect and where the cable may be stepped on.

MAINTENANCE

It does not require any maintenance other than periodic cleaning.

The laser device is intended for indoor use only.

CLEANING THE DEVICE

- Turn off the laser device and unplug from the electrical wall outlet.
- Regularly clean the laser area with a damp cloth. Do not use harsh detergents or cleaning products on any part of your device.

DISPOSAL INSTRUCTIONS

This symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment.

The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human

health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

ELECTRICAL SAFETY

Basic safety precautions should be observed when using this, or any other electrical device, especially when children are present.

- ▲ **WARNING:** Keep away from water. Contact with water may result in an electrical shock.

To reduce the risk of hazard by electric shock:

- ▲ Do not use while bathing.
- ▲ Do not place or store the laser device where it can fall or be pulled into a tub or sink.
- ▲ Do not expose the laser device to running water, or drop into water or any other liquid. If the laser device falls into water, unplug it immediately from the power outlet.
- ▲ If using the laser device near water, the use of an electrical outlet protected by a Ground Fault Circuit Interrupter is recommended.
- ▲ Do not dispose of the laser device in a fire or store in high temperatures greater than 140°F (60°C).
- ▲ Do not turn on or use the laser device while it is charging.
- ▲ If you have any problems with this device, such as setting up, maintaining or using, please contact customer service.
- ▲ Please report to customer service if any unexpected operation or events occur.
- ▲ After drop/shock, that may cause changes in the performance, please contact customer service.
- ▲ Be careful, do not put the power cord around your neck, it may cause strangulation.
- ▲ Keep the device out of the reach of children/pets to avoid inhalation or swallowing of small parts.
- ▲ Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- ▲ The device should be used only with the accessories recommended for use by the manufacturer.
- ▲ When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight.

ELECTROMAGNETIC EMISSIONS GUIDANCE**Guidance and manufacturer's declaration – electromagnetic emissions**

The Low Level Laser Therapy (LLLT) device is intended for use in the electromagnetic environment specified below.

The customer or the user of the laser device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The laser device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The laser device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliance	
<p>This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:</p> <ul style="list-style-type: none"> - Reorient or relocate the receiving antenna. - Increase the separation between the equipment and receiver. - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. - Consult the dealer or an experienced radio/TV technician for help. 		
<p>This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) the device must accept any interference received, including interference that may cause undesired operation.</p>		

ELECTROMAGNETIC IMMUNITY GUIDANCE

Guidance and manufacturer's declaration – electromagnetic immunity

The Low Level Laser Therapy (LLLLT) device is intended for use in the electromagnetic environment specified below.

The customer or the user of the laser device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines Non applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	± 1kV differential mode Non applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the laser device require continued operation during power mains interruptions, it is recommended that the laser device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The laser device power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

ELECTROMAGNETIC IMMUNITY GUIDANCE

Guidance and manufacturer's declaration – electromagnetic immunity

The Low Level Laser Therapy (LLLT) device is intended for use in the electromagnetic environment specified below.

The customer or the user of laser device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 KHz to 80 MHz 3 V/m 80MHz to 2,5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the laser device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the laser device is used exceeds the applicable RF compliance level above, the laser device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the laser device.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

ELECTROMAGNETIC DISTANCE RECOMMENDATIONS

Recommended separation distance between portable and mobile RF communications equipment and the Low Level Laser Therapy (LLLT) device.

The laser device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the laser device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the laser device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

WARRANTY

ANY ATTEMPT TO SELF-SERVICE THE LASER DEVICE, UNLESS INSTRUCTED BY HAIRMAX INTERNATIONAL CUSTOMER SERVICE, WILL VOID YOUR WARRANTY

To validate your warranty, please register online at: www.hairmax.com/warranty

THE DEVICE IS WARRANTED BY HAIRMAX INTERNATIONAL, LLC ("HAIRMAX") AGAINST DEFECTS IN MATERIALS OR WORKMANSHIP FOR 12 MONTHS FROM DATE OF PURCHASE. IF THIS PRODUCT IS DETERMINED TO BE DEFECTIVE, HAIRMAX WILL REPAIR OR REPLACE IT, AT ITS OPTION, AT NO CHARGE DURING THE WARRANTY PERIOD. THE WARRANTY PERIOD IS NOT EXTENDED IF HAIRMAX REPAIRS OR REPLACES THE PRODUCT. AFTER THE WARRANTY PERIOD, THE USER MUST PAY FOR ALL PARTS AND LABOR CHARGES. IF THE PRODUCT IS REPAIRED OR REPLACED, EITHER WITHIN OR OUTSIDE THE WARRANTY PERIOD, HAIRMAX RESERVES THE RIGHT TO PROVIDE EQUIVALENT PARTS IN DIFFERENT COLORS TO YOUR ORIGINAL PRODUCT.

THIS WARRANTY DOES NOT COVER DAMAGES CAUSED BY ACTS OF GOD, MISUSE, NEGLIGENCE, ACCIDENT, MODIFICATION OF OR TO ANY PART OF THE PRODUCT. THIS WARRANTY DOES NOT COVER DAMAGE DUE TO IMPROPER OPERATION OR MAINTENANCE, CONNECTION TO IMPROPER VOLTAGE SUPPLY, OR ATTEMPTED REPAIR BY ANYONE OTHER THAN HAIRMAX. THIS WARRANTY IS VOID IF PRODUCT IS TAMPERED WITH.

TO OBTAIN WARRANTY SERVICE, YOU MUST CALL HAIRMAX INTERNATIONAL CUSTOMER SERVICE TO REQUEST AN AUTHORIZATION NUMBER. YOU MUST DELIVER THE PRODUCT FREIGHT PREPAID, IN EITHER ITS ORIGINAL PACKAGING OR PACKAGING AFFORDING AN EQUAL DEGREE OF PROTECTION, TOGETHER WITH YOUR NAME, ADDRESS AND AUTHORIZATION NUMBER, TO:

Hairmax International, LLC, 1040 Holland Drive, Boca Raton, FL 33487, USA






HAIRMAX INTERNATIONAL, LLC SHALL NOT BE LIABLE FOR ANY SPECIAL, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES THAT ARE DIRECTLY OR INDIRECTLY RELATED TO THE USE OF, OR THE INABILITY TO USE, THE PRODUCT, EVEN IF HAIRMAX INTERNATIONAL, LLC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES WHETHER IN AN ACTION OF CONTRACT, NEGLIGENCE OR OTHER TORTUOUS ACTION. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF SPECIAL, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU. THIS WARRANTY IS THE ONLY WRITTEN OR EXPRESS WARRANTY GIVEN BY HAIRMAX INTERNATIONAL, LLC. ALL IMPLIED WARRANTIES FOR THE PRODUCT, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN TIME TO THE TERM OF THIS WARRANTY. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU.

GOVERNING LAW; JURISDICTION AND VENUE: THIS WARRANTY, AND THE TERMS OF SALE OF THE PRODUCT, SHALL BE GOVERNED BY AND CONSTRUED IN ACCORDANCE WITH THE LAWS OF THE STATE OF FLORIDA WITHOUT REGARD TO PRINCIPLES OF CONFLICTS OF LAWS. THE PARTIES IRREVOCABLY SUBMIT AND CONSENT TO THE SOLE AND EXCLUSIVE JURISDICTION AND VENUE OF THE FLORIDA STATE COURTS IN AND FOR PALM BEACH COUNTY, FLORIDA, U.S.A. AND THE FEDERAL COURTS IN AND FOR THE SOUTHERN DISTRICT OF FLORIDA, U.S.A. THE PARTIES AGREE NOT TO RAISE THE DEFENSE OF FORUM NON CONVENIENS.

WARRANTY OUTSIDE THE U.S.A.

THE WARRANTY ON THIS PRODUCT SHALL COMPLY WITH APPLICABLE LAW WHEN SOLD BY AN AUTHORIZED DEALER TO A CONSUMER OUTSIDE OF THE UNITED STATES, INCLUDING EXTENDED WARRANTY COVERAGE PERIODS AS MAY BE LEGALLY REQUIRED IN SUCH JURISDICTIONS.

SYMBOLS

	Manufacturer		FCC compliance mark
	Date of manufacture		RoHS
	Country of manufacture		Laser warning
	Serial #		Low-level ESD
	Keep Dry		Caution laser radiation class 3R
	Temperature limit		Stand-by
	Humidity limitation		On/Off (push-push)
	Caution		Direct current
	CE mark European conformity		EC Rep
	General warning sign		AU Rep
	Refer to instruction manual/ booklet		BR Rep
	WEEE		UK Rep
	Atmospheric pressure limitation		Laser on
	RF transmitter		Laser beam; radiation of laser apparatus
IP22 Ingress Protection Code			For indoor use only
	Type BF applied part		InMetro
†LASER APERTURE† Laser aperture			Société Générale de Surveillance

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Copyright 2025 Hairmax Intl., LLC. All Rights Reserved.
Emergo Consulting (UK) Limited c/o