bridge

Vibroacoustic Massage Cushion

Recharge Body. Elevate Mind.

Instruction Manual

BVA-HB150-EB

Experience vibroacoustic therapy for both body and mind with

bridge

Sending low frequency soundwaves deep into the body, relaxing tight muscles, calming the nervous system and restoring mind-body balance.

SET UP:

- Slide the correct pin attachment (EU 2-pin or GB 3-pin) onto the adapter plug until it clicks into place.

- Attach the adapter to the cushion by connecting the jack socket, then plug the adapter into a wall socket.
- Place the cushion either flat on a bed or sofa to lay on it, or against the back of a comfy chair to use it seated. We recommend that you also dim the lights and ensure you are undisturbed to get the best possible experience.

USING THE REMOTE CONTROL:

- Press the power button  to turn on.
- Press the mode button  to select a track.
- Press the timer button  to change the duration: 30 min, 20min, 10min or 5 min
(Note: The times are for guidance only as all tracks vary slightly in length)
- Use the + and - buttons to adjust vibration intensity () and volume ()
- Press the mute button  to reduce the sound to its lowest level. (Please just note that because the unit is powered by sound, it will not mute completely, or the vibration would stop too)

STREAMING MUSIC VIA BLUETOOTH:

Explore other vibroacoustic tracks by streaming them from your phone/tablet to **bridge**.

- Press the Bluetooth button  on the remote control to put the product into 'connect' mode (white light flashing)
- On your phone (or tablet), in Bluetooth settings, search for and connect to **bridge vibroacoustic cushion**
- Search for vibroacoustic tracks in Spotify / Apple Music etc and play them through the product, just as you would with an external speaker.
- We have also crafted a bespoke music list for you which you can find by searching for **bridge vibroacoustic** by Homedics on Spotify, Apple Music or YouTube Music

AUTO SWITCH OFF:

If the product is not connected to Bluetooth and has been left idle for 15 mins, it will automatically switch off.

CLEANING AND CARE:

- Always unplug the product before cleaning.
- Clean the surface of the product with a damp sponge only. Do not saturate. Allow to dry fully before use.
- Avoid the use of chemical or abrasive cleaners which can damage the product.

**READ ALL INSTRUCTIONS BEFORE USE.
SAVE THESE INSTRUCTIONS FOR FUTURE
REFERENCE.**

- ALWAYS unplug the appliance from the electrical outlet immediately after using and before cleaning. To disconnect, turn all controls to the 'OFF' position, then remove plug from outlet.
- NEVER leave an appliance unattended when plugged in. Unplug from the mains outlet when not in use and before putting on or taking off parts or attachments.
- This appliance is suitable for ages 8 and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Do not place or store appliance where it can fall or be pulled into a bath or sink. Do not place in or drop into water or other liquid.
- DO NOT reach for an appliance that has fallen in water or other liquids. Switch off at the mains and unplug it immediately. Keep dry – DO NOT operate in wet or moist conditions.
- NEVER insert pins, metallic fasteners or objects into the appliance or any opening.
- Use this appliance for the intended use as described in this booklet. DO NOT use attachments not recommended by Homedics.
- NEVER operate the appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return it to the Homedics Service Centre for examination and repair.
- Keep cord away from heated surfaces.
- DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
- DO NOT operate under a blanket or pillow. Excessive heating can occur and cause fire, electrocution or injury to persons.
- DO NOT carry this appliance by the cord or use cord as a handle.
- DO NOT use outdoors.
- DO NOT crush. Avoid sharp folds.
- This product requires a 100-240 V AC power supply.
- The appliance is only to be used with the power supply unit provided with the appliance.
- DO NOT attempt to repair the appliance. There are no user serviceable parts. All servicing of this appliance must be performed at an authorised Homedics Service Centre.
- If you have any concerns regarding your health, consult a doctor before using this appliance.
- This is a non-professional appliance, designed for personal use. DO NOT use as a substitute for medical attention.
- Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your GP.
- Not recommended for use by individuals with sensory deficiencies including diabetic neuropathy.
- DO NOT use on an infant, invalid or on a sleeping or unconscious person.
- DO NOT use on insensitive skin or on a person with poor blood circulation.

- This appliance should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls.
- NEVER use directly on swollen or inflamed areas or skin eruptions.
- Do not use during pregnancy
- DO NOT use if you have a pacemaker or implanted defibrillator
- DO NOT use if you suffer from severe osteoporosis, thrombosis, varicose veins, unstable blood pressure, heart conditions or risk of blood clots.
- Avoid using the product on areas of acute inflammation, infection, or recent fractures, surgery or injury.
- Failure to follow the above may result in the risk of fire or injury.

3 YEAR GUARANTEE

FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 3 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK / EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and / or authorised, or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages.

To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (as proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through Homedics Service Centre. Service of this product by anyone other than the Homedics Service Centre voids the guarantee. This guarantee does not affect your statutory rights.

For your local Homedics Service Centre, go to

www.homedics.co.uk/servicecentres

WEEE explanation

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

Operating frequency bands: 2402 - 2480MHz

CE BT Max Power Transmitted Power: 1.60 dBm

In off mode, this product will consume power at a rate of 0.03W.

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FKA Brands Ltd, Somerhill Business Park, Tonbridge, Kent TN11 0GP, UK
FKA Brands Ltd, 29 Earlsfort Terrace, Dublin D02AY28, Ireland