# FACELIFT & TONE

## BEAUTY SPA

 the beauty specialists

ENGLISH	4
FRANÇAIS	13
DEUTSCH	22
NEDERLANDS	31
ESPAÑOL	40
PORTUGUÊS	49
ITALIANO	58

## WARNINGS

- **Please read all instructions thoroughly and retain for future reference.**
- This appliance can be used by children aged 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and maintenance shall not be made by children without supervision.
- Do not place or store in direct sunlight. Store in a cool dry place.
- Keep out of the reach of children and pets.
- If you suffer from a skin condition, consult your dermatologist before use
- For hygiene reasons we recommend for single person use only.
- A safe charging temperature range would be 0-40°C.
- Do not use this appliance on loose skinned parts of your body or any areas showing signs of: varicose veins, skin diseases (for example eczema, psoriasis, dermatitis), cuts and abrasions, irritation, warts, moles, skin tags, skin rash, hypersensitivity or sunburn, phlebitis, recent scar tissue, other skin conditions. We recommend you consult your physician before treatment if you suffer from any of the above conditions.
- Do not use on cuts or broken skin.
- Do not use on areas that have had cosmetic surgery.
- Do not leave the device unattended. Turn off before you leave the room.
- This appliance must only be supplied at safety extra low voltage (SELV) corresponding to the marking on the appliance.
- Do NOT use this product if it is not working properly, or if it has been dropped or damaged.
- Do not charge for longer than 24 hours.
- Do not use immediately around the eyes.
- Do not use on areas that have had cosmetic surgery.
- Do NOT immerse this product in water.
- Do NOT hold under running water.
- The appliance may be cleaned by wiping with a damp cloth. Do NOT use solvents or abrasive cleaners.
- For external use only.

## KNOW YOUR UNIT

(see page 3)

- |                                                    |                                              |                              |
|----------------------------------------------------|----------------------------------------------|------------------------------|
| <b>1</b> Treatment Head                            | <b>5</b> On/Off/Mode                         | <b>10</b> Conductive Contact |
| <b>2</b> Intensity/Temperature Indicator           | <b>6</b> Treatment Lights (Red, Blue, Green) | <b>11</b> USB-C Socket       |
| <b>3</b> Mode Indicator (Clean, EMS, Pulse & Cool) | <b>7</b> Magnetic Fixing Ring                | <b>12</b> Stand              |
| <b>4</b> Intensity/Temperature                     | <b>8</b> Battery Indicator                   | <b>13</b> USB Charging Cable |
|                                                    | <b>9</b> LED Light                           | <b>14</b> Conductive Gel     |

## INTRODUCTION

Rio FaceLift & Tone Beauty Spa combines six technologies inspired by professional treatments - warming, cooling, EMS microcurrent, EMS intermittent pulses, sonic vibration, and tri-colour LED Light therapy - each offering targeted benefits for the skin.

Combined into one time-efficient, holistic device, these technologies deliver a complete, step-by-step, spa-inspired treatment to deeply cleanse, lift, enhance product absorption, and lock in hydration.

Whether you're looking to soften the visible signs of ageing, smooth skin texture, or maintain a healthy, radiant glow, Rio FaceLift & Tone Beauty Spa is designed to follow the natural contours of the face, delivering a complete, results-oriented skincare protocol in the comfort of your own home.

## THE SIX TECHNOLOGIES

**WARMING:** Helps relax facial tension and soften the skin's surface to enhance cleansing and product absorption. The heat function (from 38°C up to 42°C) supports collagen production and cellular repair while gently opening pores for a deeper cleanse.

A Thermal Control System ensures a safe, effective temperature - protecting skin and preserving fat cells.

**SONIC VIBRATION:** Operating at 8,000 vibrations per minute, stimulates bloodflow and helps activate lymphatic drainage, encouraging the natural elimination of toxins and promoting overall skin vitality.

**EMS - MICROCURRENT:** Low-level electrical impulses gently stimulate the facial muscles, supporting toning, firming and lifting. Muscles in the face are unique in that they are the only muscles in the body which connect to the skin, making EMS facial toning particularly effective.

**EMS - INTERMITTENT PULSES:** Delivers rhythmic pulses inspired by traditional acupressure techniques to help relieve muscle tension and soften the appearance of lines and tired skin. Stimulates microcirculation, bringing oxygen and nutrients to the skin, to help revitalise dull-looking skin.

**COOLING:** Provides an instant soothing effect (from 15°C down to 6°C ), helping to reduce puffiness, refine pores, and calm stressed skin. Locks in moisture, enhances firmness, and leaves the skin feeling fresh, smooth, and toned.

**LED LIGHTS:** **Blue Light:** Targets the surface and near surface, to promote clear skin and an even complexion. **Green Light:** Helps promote skin recovery and enhance natural radiance. **Red Light:** Penetrates further into the skin to support cellular metabolism and collagen production.

## IMPORTANT - PATCH TEST

Before you start you must ALWAYS do a patch test with the Conductive Gel on an inconspicuous area of your skin for possible skin reaction. Wait 24 hours. If irritation occurs do not continue treatment.

**Note:** Always start with the lowest level for each mode to allow your skin to adjust. You can gradually increase the intensity as you become more comfortable.

The sonic vibration function activates on skin contact when your hand is also touching the Conductive Contacts (10) on the sides of the unit. It stops when released, saving energy and extending battery life.

## CHARGING

- Connect the unit to a USB port using the USB Charging Cable.
- The Battery Indicator will flash while the device is charging. A single bar indicates a low battery level, while 3 bars signify a full charge. Once fully charged, the device will automatically power off after 5 minutes.
- It takes about 3 hours to fully charge.
- When fully charged the device should operate for approximately 40 minutes.

**IMPORTANT:** Disconnect from USB Charging Cable before using.

## USING YOUR FACELIFT & TONE

1. Clean and dry the skin in the area to be treated.
2. Press the On/Off/Mode button for 2 seconds to turn on the FaceLift & Tone.  
The default "CLEAN" mode will activate automatically.
3. To select a different treatment mode, press the Mode button to cycle through the four options:

**CLEAN:** Deep, relaxing cleanse

**EMS:** Lift and tone your skin

**PULSE:** Boost product absorption

**COOL:** Lock in moisture, refresh skin

4. Depending on the treatment mode, use the Intensity/Temperature button to select the following:

**CLEAN:** Warming temperature: 38°C, 40°C, 42°C

**EMS:** Microcurrent intensity: Low, Medium, High

**PULSE:** Pulses intensity: Low, Medium, High

**COOL:** Cooling temperature: 15°C, 10°C, 6°C

The selected intensity/temperature level is indicated (2).

5. To choose the LED Light colour, press the LED Light button to cycle through the available options: No light, Red, Blue, or Green. LED selection is available in every mode.
6. Move the Treatment Head slowly over the skin, following the direction of the arrows (see Treatment Areas). Each movement should last 5 to 8 seconds. Repeat 4 to 5 times per area.
7. To turn off the device, press and hold the On/Off/Mode button for 2 seconds.

### Important:

- Do not use for more than 5 minutes in any Treatment Area.
- The unit will automatically turn off after 5 minutes for the CLEAN, EMS and PULSE treatment modes, and after 2 minutes for the COOL mode.

# TREATMENT AREAS

① Forehead    ② Right Face    ③ Left Face    ④ Neck

## DISCOVER THE 4 COMBINED TREATMENTS

Each treatment can be carried out separately or in combination in following sequence:

### **CLEAN** Warming + Sonic Vibration + Light Therapy

- *The best start to your beauty routine*

Gently detach the Magnetic Fixing Ring from the device. Position a cotton pad over the Treatment Head. Reattach the ring to hold it securely.

Apply your favourite cleansing water or makeup remover. Glide the tool slowly from the centre of the face outwards to deeply cleanse and prepare the skin for the next mode.

**Tip:** If you enjoy double cleansing, try this protocol. Start with an oil-based product to dissolve makeup, sunscreen, sebum, and other oil-based impurities. Follow with a water-based cleanser.

### **EMS** Microcurrent + Sonic Vibration + Light Therapy

Apply a generous amount of Conductive Gel to the treatment head or directly on to the skin, one area at a time. Glide the tool slowly from the centre of the face outwards. If the gel starts to dry, reapply.

**Tip:** EMS is more effective when your skin is well hydrated—drink a glass of water before the session and/or apply a water-based hydrating serum before the gel. You can also use Rio Hydrogel Mask (available at [www.riobeauty.co.uk](http://www.riobeauty.co.uk)) instead of the Conductive Gel.

### **PULSE** Intermittent pulses + Sonic Vibration + Light Therapy

- *Can be used post-cleansing and/or after EMS treatment.*

Apply your favourite serum, essence, eye cream or use over a tissue mask to boost product absorption.

**Tip:** Move slowly and gently across the face - this allows the micro-pulses to maximise absorption without dragging the product.

### **COOL** Cooling + Sonic Vibration + Light Therapy

- *Ideal as the final touch of your routine.*

Apply after product absorption, using a gentle press-and-hold motion, moving area by area: rest the tool on one area for a few seconds, then move slowly from the centre of the face outwards.

**Tip:** Diffuses a fine mist of thermal water to enhance the refreshing sensation and comfort of the treatment.

**You can customise each treatment by selecting the LED Light colour that best suits your needs.**

## AFTER USE

- If you have very sensitive skin, you may notice redness after use. This should subside after 30 minutes.

## CLEANING AND MAINTENANCE

- Before and after use, wipe the Treatment Head with a clean damp cloth, do not immerse the unit in water.

## DISPOSAL

The Waste of Electrical and Electronic Equipment (WEEE) Directive (2012/19/EU) has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill. This product is classed as Electrical or Electronic equipment so please ensure that at the end of the life of this product it is disposed of in the correct manner in accordance with local authority requirements. It must NOT be disposed of with household waste.

### **Recycle Responsibly:**

1. Check your local municipal guidelines
2. Separate the components and dispose of them correctly

## SPECIFICATION

Model: FCTL  
Input: 5.0V  $\equiv$  1.0A 5.0W Max   
Battery: Li-Po 3.7V 1000mAh 3.7Wh

Read these instructions carefully before use and then retain for future use.

This product is for indoor use only.

This product complies with current EU product legislation.

This unit complies with current UKCA product legislation.

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© The Dezac Group Ltd 2025 Manufactured by The Dezac Group Ltd PO Box 17, Cheltenham, England GL53 7ET