hairmax®
TRANSFORMATIVE HAIR REGROWTH

QUICK START GUIDE

LASERBAND® 82

Hair Growth **Laser Light Device**

ComfortFlex

QUICK START GUIDE

LaserBand 82 *comfortflex*

3 times
per week


FOR MAXIMUM BENEFIT

- :: Scalp should be freshly washed (hair can be dry or damp).
- :: Use the laser device prior to applying styling products (these can build up and impede the laser light).
- :: Regularly clean the laser area with a damp cloth. Do not use harsh detergents or cleaning products on any part of your device.


CHARGING OPTIONS

You must charge the device for 4-6 hours before initial use


Prior to charging the laser device, make sure the device is turned off.

OPTION 1 Plug the power cord into the charging stand and the other end into an electrical wall outlet. Place the device on the charging stand. The device indicator white light will flash during the charging process. A steady light indicates the device is fully charged.

OPTION 2 Place the power cord directly into the charging socket on the side of the device and the other end into an electrical wall outlet.

The laser device indicator white light will flash during the charging process. A steady light indicates the device is fully charged.

Do not use the laser device while charging.

Prior to treatment, test your skin for sensitivity by shining the light on the inside of your forearm for 2 minutes. Make sure your forearm is clean before testing. If you see any reaction, do not use your laser device and contact customer service.

Before Starting Your Treatment:

1. Make sure the battery is FULLY CHARGED.
2. Test your Hairmax® laser device for sensitivity.

TREATMENT METHOD

The LaserBand will be used in 3 positions on your scalp. Leave it in place for a minimum of 30 seconds for each position. It will vibrate every 30 seconds to let you know when to move it to the next position.

To start the treatment, gently press the power button once to activate the red lasers.

1. Place the laser device at your hairline and leave it in place until it vibrates. This is position 1.


2. When you feel the vibration, slide the laser device backwards about 3 inches (7.5 cm) to the middle portion of your scalp.

Use light pressure so the teeth can separate your hair. Leave it in position 2 until you feel the vibration again.

3. Repeat the 30 second procedure on position 3. You have now completed today's treatment.

After treatment, remove the laser device from your head and turn off the device by pressing the power button once. Place the device in the charging cradle/stand and store for future use.

NOTE: The LaserBand 82 will turn off automatically after 10 minutes.

Based on the durability of the lasers and electronic components, the expected life of the laser device is designed to last up to 10 years of regular use.

For detailed information about troubleshooting, proper use, warnings and cautions, access the full user manual by scanning the QR code located on the back panel. Please read the entire user manual before operating.

The Hairmax® laser device has a limited warranty against defects in materials or workmanship from Hairmax International, LLC. Items that are damaged or abused are not covered under the warranty and will be charged at their regular repair rates. Any attempt to self-service the device, unless instructed by HAIRMAX INTERNATIONAL CUSTOMER SERVICE, will void your warranty. Please refer to our website for complete warranty terms.


**To validate your warranty, register online at
www.hairmax.com/warranty**

Laser Light - Avoid Direct Eye Exposure

"Hairmax® laser devices are indicated to promote hair growth in males with Androgenetic alopecia who have Norwood-Hamilton Classifications of IIa - V, and in females with Androgenetic alopecia who have Ludwig (Savin) Classifications I - II or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I - IV.

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SCAN QR CODE

To access the full
user manual

Or visit: www.hairmax.com/usermanual

You can also call customer service free of charge to request a hard copy

CUSTOMER SERVICE

Direct line +1 561-417-0200

United States and Canada (toll-free) 1-866-527-3726

Authorized Representative in Europe (Regulatory Affairs Only)

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Tel: +31 (0) 70 345-8570 Fax: +31 70 346-7299

Brazil Sponsor

EMERGO BRAZIL
Avenida Francisco Matarazzo 1.752,
Salas 502/503
Água Branca
São Paulo - SP CEP – 05001-200

Tel: +55 11 4800 9226

Australia Sponsor

EMERGO AUSTRALIA
201 Sussex Street
Darling Park, Tower II, Level 20
Sydney NSW 2000

Tel: +61 2 9006 1662

United Kingdom Sponsor

EMERGO UNITED KINGDOM
Emergo Consulting (UK) Limited c/o
Cr360 – UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ

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1040 Holland Drive | Boca Raton, Florida | 33487 | USA

www.hairmax.com

Copyright 2025 Hairmax Intl., LLC. All Rights Reserved.

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HAIRMAX INTERNATIONAL, LLC