
LUCENT

SPARK

Compact Electric Toothbrush

MANUAL AND WARRANTY

SPARK

USER

MANUAL

INSTRUCTIONS OF USE	ENGLISH	4-5
HANDBUCH	DEUTSCH	6-7
MODE D'EMPLOI	FRANÇAIS	8-9
MANUAL DE USUARIO	ESPAÑOL	10-11
MANUALE D'USO	ITALIANO	12-13
GEbruiksaanwijzing	NEDERLANDS	14-15
BRUKSANVISNING	SVENSKA	16-17
BRUGSANVISNING	NORSK	18-19
BRUGSANVISNING	DANSK	20-21
KÄYTTÖOHJEET	SUOMALAINEN	22-23
INSTRUKCJA UŻYCIA	POLSKI	24-25

THANK YOU FOR BUYING THE SPARK TOOTHBRUSH

Spark is the compact sonic toothbrush built for simplicity and performance. With powerful vibrations, a sleek design, and USB Type C charging, it's made to fit your life. At home or on the go.

HOW TO USE SPARK

On/Off:

Press the button briefly to turn the toothbrush on or off.

Lock/Unlock:

Press and hold the button for 5 seconds. When locking, all three indicator lights flash once. When unlocking, all three lights blink three times.

Change Mode:

While the toothbrush is on, press and hold the button for 2 seconds to switch modes.

Replace Brush Head:

Gently pull the old brush head off the handle. Press a new brush head firmly onto the metal pin. For the best cleaning results, replace the brush head regularly.

Timer Function:

The toothbrush pauses briefly every 30 seconds and stops after 2 minutes to help guide effective brushing time.

Charging:

The USB Type C connection is placed under the lid at the bottom of the toothbrush. Put the other side of the cable into a USB Type C power connection, e.g. a phone charger adapter or your laptop output. Always keep the toothbrush off when charging. The battery is fully charged after 3 hours. It's recommended to fully charge the first time or has not been used for a few months.

MODES

Sensitive:

Gentle vibrations designed to protect gums—ideal for sensitive teeth and daily use.

Power:

More intense sonic vibrations for a deeper, more thorough clean.

Whiten:

Pulsating vibrations help remove plaque and surface stains for a brighter smile.

PRECAUTION

- Keep out of reach of children under 3 years of age.
- Do not use if the toothbrush or charger appears damaged.
- Only use with a compatible power source (max 5V input).
- Do not expose to temperatures above +60°C or below -20°C.
- Replace the brush head at least every 3 months for optimal hygiene.
- Do not submerge the handle or charger in water.
- Avoid using while charging.

PRODUCT GUIDE AND SPECS

- | | |
|------------------------|--------------------|
| 1. Soft silk bristles | |
| 2. Vibration frequency | 35,400 RPM |
| 3. One button control | |
| 4. Mode indicators | |
| 5. Water proof | (IPX7) |
| 6. USB-C connection | |
| 7. Li-on battery: | 3.7V 330mAh |
| 8. Working power: | 1.8W |
| 9. Input voltage: | 5V – 1A |
| 10. Charging time: | Ca 3 hours |
| 11. Battery life: | Ca 3 weeks |
| 12. Charging input: | Max 5V |

WARRANTY

Our products are made to meet high standards and deliver a great experience. If a manufacturing defect occurs, we offer a one-year warranty from the date of purchase covering defects in materials and workmanship.

To claim the warranty, please present the original purchase receipt, to either the point of purchase or to our customer service.

The warranty does not cover misuse, normal wear and tear, unauthorised repairs or damage caused by external events.

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Be Lucent AB