

OPERATOR: Helena Lamplough	DATE: 21.07.2024 09:21	PROCESS C/WN	
CUSTOMER: SBC / Tricogen - SBC	ORDER NUMBER: BM06958	PROCESS MAGENTA	
UNIQUE NUMBER: LSBCS150	ITEM CODE: GRG425PB	PROCESS YELLOW	
PRINTER: All4Labels Leeds		PROCESS BLACK	
CUTTER SSL: x120x165mm_rectangle_3mm_rads_300ml_wrap-02		PANTONE 369 C	
CUTTER IN BRANDWORKS:		UV WHITE1	
CUTTER NOTES: 120mm x 165mm		CUTTER	
BARCODE DATA: 5034127425205			
FRONT QR CODE:			
BACK QR CODE:			
CODE 1 TYPE: EAN-13			
CODE 1 BWR: 0.06 mm	CODE 1 MAG: 80%		
CODE 1 DATA: 5034127425205			
CODE 2 TYPE:			
CODE 2 BWR:	CODE 2 MAG:		
CODE 2 DATA:			
CODE 3 TYPE:			
CODE 3 BWR:	CODE 3 MAG:		
CODE 3 DATA:			
CODE 4 TYPE:			
CODE 4 BWR:	CODE 4 MAG:		
CODE 4 DATA:			
EYEMARK / CAP COLOUR:			
CHECKED BY:			

120 mm

165 mm

Hydrate · Energise · Glow

SBC Skincare's signature lightweight gel formula

A revitalising moisturising gel infused with Grapefruit and Ginger extracts to help brighten and smooth the skin on your face and body. The antioxidant-rich blend of Fennel, White Nettle and Elderflower soothe, balance and help even skin tone. Ideal for dull, dehydrated or combination skin.

Directions: Apply to clean skin as a versatile moisturiser, primer or mask. The gel can be used with electrical beauty tools on both positive (+) and negative (-) settings (always follow manufacturers' instructions).

(IT) Pompelmo e Zenzero Gel Idratante

Idrata · Energizza · Illumina

Istruzioni per l'uso: Applicare sulla pelle pulita come versatile idratante, primer o maschera. Il gel può essere utilizzato con strumenti di bellezza elettrici con impostazioni sia positive (+) che negative (-) (seguire sempre le istruzioni del produttore).

(DE) Grapefruit und Ingwer Hautpflege Gel

Befeuchtend · Erwecken · Strahlen

Anwendung: Als vielseitige Feuchtigkeitscreme, Grundierung oder Maske auf die gereinigte Haut auftragen. Das Gel kann mit elektrischen Beauty-Geräten sowohl in positiver (+) als auch in negativer (-) Einstellung verwendet werden (beachten Sie stets die Anleitungen des Herstellers).

Ingredients: Aqua, Alcohol Denat., Glycerin, Citrus Grandis (Grapefruit) Peel Oil, Zingiber Officinale (Ginger) Extract, Lamium Album (White Nettle) Flower/Leaf/Stem Extract, Foeniculum Vulgare (Fennel) Seed Extract, Sambucus Nigra (Elderflower) Flower Extract, Triethanolamine, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Cocamidopropyl Betaine, Polysorbate 20, Ethylhexylglycerin, Tetramethyl Acetyloctahydronaphthalenes, Sodium Chloride, Phenoxyethanol, Parfum, Limonene, Linalool, Citrus Aurantium (Orange) Peel Oil, Hexamethylindanopyran, Terpineol, Geranyl Acetate, Geraniol, Vanillin, Citral, Citric Acid, Sodium Benzoate, Citronellol, Potassium Sorbate, CI 16035, CI 15985, CI 42090.

Due to natural ingredients, colour may vary or change over time.

1

PET

BOTTLE

BOTTIGLIA

FUSIONE

Widely Recycled

Collezione di Plastica

Plastiksammlung

92

C/PP

PUMP

POMPA

POMPE

Not Currently Recycled

Non Ordinato

Unsortiert

MADE IN THE UK

SBC Skincare LTD RH1 2NB Redhill UK

www.sbcskincare.co.uk

Cura Marketing GmbH 6020 Innsbruck Austria

www.curacosmeticsgroup.com

GRG425PB/01

12M

5 034127 425205 >

British Brand

In The UK

30% PCR

Recycled Plastic

Vegan Friendly

SBC

SKINCARE

Since 1990

Grapefruit & Ginger

Moisturising Gel

Ionized + & -

Face & Body

Uplift & Energise

300ml e 10.1 US Fl. Oz.

SSL REF:
120x165mm_rectangle_3mm_rads_300ml_wrap - 02

OPERATOR: Helena Lamplough	DATE: 09.07.2024 11:44
CUSTOMER: SBC / Tricogen - SBC	ORDER NUMBER: BM06827
UNIQUE NUMBER: TTRIC075	ITEM CODE: GRG432FT-01
PRINTER: Mpack Poland Sp. z o.o.	
CUTTER SSL: Mpack_tube_FI40x120-02	
CUTTER IN BRANDWORKS:	
CUTTER NOTES: F140X120	
BARCODE DATA: 5034127432111	
FRONT QR CODE:	
BACK QR CODE:	
CODE 1 TYPE: EAN-13 CODE 1 BWR: 0 mm CODE 1 DATA: 5034127432111	CODE 1 MAG: 80%
CODE 2 TYPE: CODE 2 BWR: CODE 2 DATA:	CODE 2 MAG:
CODE 3 TYPE: CODE 3 BWR: CODE 3 DATA:	CODE 3 MAG:
CODE 4 TYPE: CODE 4 BWR: CODE 4 DATA:	CODE 4 MAG:
EYEMARK / CAP COLOUR:	
CHECKED BY:	



varnish area, varnish to be added by the printer

SSL reference: Mpack_tube_FI40x120-02